# Comparison of EVAR and ESAR for Infrarenal Aortic Aneurysms With a Wide Proximal Neck

Status: RECRUITING

## Eligibility Criteria

Age: 18 years and over

This study is NOT accepting healthy

Healthy Volunteers: volunteers

#### **Inclusion Criteria:**

\* Age 18 years or older \* Provided written informed consent \* Clinical necessity for treatment of the AAA, according to the current guidelines in the geographies participating \* Elective repair \* Eligible anatomy for treatment with the Endurant II/IIs stent graft system and Heli-FX EndoAnchor system according to the IFU of both devices \* Infrarenal neck diameter ≥ 28 mm and ≤32 mm \* Proximal neck length ≥10mm

### **Exclusion Criteria:**

\* Anatomy outside the IFU of the ndurant II/IIs stent graft system and Heli-FX EndoAnchor system \* Planned use of AUI main body device \* Patient is participating in another clinical study, potentially conflicting with the outcomes of the current study. \* Patient with eGFR \< 30 ml/min/1.73m2 before the intervention \* Patient's life expectancy \<2 years as judged by the investigator \* Patient has a psychiatric or other condition that may interfere with the study \* Patient has a known allergy to any device component \* Patients with a systemic infection who may be at increased risk of endovascular graft infection. \* Patient has a coagulopathy or uncontrolled bleeding disorder \* Patient has a ruptured, leaking, or mycotic aneurysm \* Patient is not eligible for standard EVAR \* Patient had a Cerebro Vascular Accident (CVA) or a myocardial infarction (MI) within the prior three months \* Patient is pregnant (Female patients of childbearing potential only) \* Patient has active COVID-19 infection or has been diagnosed with long COVID-19 requiring hospitalization within the 6 months prior to procedure. \* Patient has previously been treated with stent grafts in the aorto-iliac arteries

## Conditions & Interventions

### Interventions:

DEVICE: endovascular aneurysm repair (EVAR) using the Endurant II/IIs stent graft system, DEVICE: endosuture aneurysm repair (ESAR) using Heli-FX EndoAnchor system

#### Conditions:

Aortic Aneurysm, Abdominal

### Keywords:

AAA, Endovascular treatment, Wide proximal neck, EndoAnchors, Infrarenal, EVAR (Endovascular aneurysm repair), ESAR (Endosuture aneurysm repair), Endurant II/IIs stent graft system, Heli-FX EndoAnchor system

### More Information

Contact(s): Suzanne Holewijn, PhD - herculestrial@rijnstate.nl

**Principal Investigator:** 

Phase: NA IRB Number:

**System ID:** NCT05484115

Thank you for choosing StudyFinder. Please visit http://studyfinder.cctr.vcu.edu to find a Study which is right for you and contact ctrrecruit@vcu.edu if you have questions or need assistance.